


Anti-Cd271 Rabbit Polyclonal Antibody For Ihc-P And Elisa - P08138

Item Number: CM0020503



Introduction

High-quality rabbit polyclonal antibody targeting human CD271/p75NTR for IHC-P and ELISA applications. This antibody detects the low-affinity neurotrophin receptor p75NTR, a key regulator of neuronal development, apoptosis, and circadian rhythm. Ideal for neuroscience and apoptosis research.

[Learn More](#)

Parameter	Specification
Product Name	CD271 Rabbit pAb
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR; p75NTR/NGFR
Species	Human
GeneID (Human)	4804
GeneID	4804
Immunogen	Recombinant protein corresponding to a sequence within amino acids 330-427 of human CD271 (NP_002498.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08138
Protein Weight	45kDa
Shipping	Ice bag