

Anti-Cd272/Btla Monoclonal Antibody For Fc - Q7Tsa3

Item Number: CM0008732



Introduction

Rabbit monoclonal antibody against mouse CD272/BTLA, conjugated to . Validated for flow cytometry (FC) to analyze CD272/BTLA expression in mouse lymphocytes. Research use for inhibitory receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD272/BTLA mAb
Remarks/Alias	A630002H24
Species	Mouse
GeneID (Human)	151888
GeneID	208154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-176 of mouse CD272/BTLA (NP_001032808.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q7TSA3
Protein Weight	34kDa
Shipping	Ice bag