

## Anti-Cd272/Btla Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q7Tsa3

Item Number: CM0006782



### Introduction

Recombinant rabbit monoclonal antibody against mouse CD272/BTLA, an inhibitory receptor on lymphocytes. Validated for WB, IF/ICC, FC, and ELISA. Ideal for studying immune checkpoint signaling.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD272/BTLA Rabbit mAb                                                                                                     |
| Remarks/Alias    | A630002H24; CD272/BTLA                                                                                                    |
| Species          | Mouse                                                                                                                     |
| GeneID (Human)   | 151888                                                                                                                    |
| GeneID           | 208154                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 30-176 of mouse CD272/BTLA (NP_001032808.2) |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | WB, IF/ICC, FC, ELISA                                                                                                     |
| Cross Reactivity | Mouse                                                                                                                     |
| SWISS            | Q7TSA3                                                                                                                    |
| Protein Weight   | 34kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |