


Anti-Cd274 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ep73

Item Number: CM0020182



Introduction

High-quality anti-mouse CD274/PD-L1 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for immune checkpoint and tumor immunology research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CD274 Rabbit pAb
Remarks/Alias	B7h1; Pd1l; Pdcd1l1; Pdcd1lg1; A530045L16Rik; 74
Species	Mouse
GeneID (human)	29126
GeneID	60533
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9EP73
Protein Weight	33kDa
Shipping	Ice bag