

Anti-Cd275 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Jhj8

Item Number: CM0019528

Introduction

Rabbit polyclonal antibody targeting human CD275 (ICOS ligand), cross-reactive with mouse. Validated for WB and ELISA. Supports research in immune regulation, T-cell biology, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD275 Rabbit pAb
Remarks/Alias	B7h; GI50; GL50; B7-H2; LICOS; B7RP-1; GL50-B; ICOS-L; Icoslg; Ly115l; mKIAA0653; CD275
Species	Human
GeneID (Human)	50723
GeneID	50723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-270 of human CD275. (NP_056605.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9JHJ8
Protein Weight	36kDa
Shipping	Ice bag