

Anti-Cd275 Rabbit Polyclonal Antibody For Wb - 075144

Item Number: CM0020879



Introduction

Anti-CD275 Rabbit Polyclonal Antibody (ICOSL, B7-H2) validated for Western blot. Detects human and mouse ICOS ligand, with a predicted molecular weight of ~33 kDa. Ideal for studies on T-cell costimulation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	CD275 Rabbit pAb
Remarks/Alias	B7h; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
Species	Human
GeneID (Human)	23308
GeneID	23308
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 19-256 of human ICOSL (NP_056074.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Mouse
SWISS	O75144
Protein Weight	33kDa
Shipping	Ice bag