

Anti-Cd276 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5Zpr3-2

Item Number: CM0023196

Introduction

High-quality rabbit polyclonal antibody against human CD276 (B7-H3), validated for Western blot, IHC-P, and ELISA. Ideal for immune checkpoint and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	CD276 Rabbit pAb
Remarks/Alias	B7H3; B7-H3; B7RP-2; 4Ig-B7-H3; CD276
Species	Human
GeneID (Human)	80381
GeneID	80381
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 29-245 of human CD276 (NP_079516.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q5ZPR3-2
Protein Weight	33kDa
Shipping	Ice bag