

## Anti-Cd276/B7-H3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5Zpr3

Item Number: CM0021141

### Introduction

Rabbit polyclonal antibody against human CD276/B7-H3 (B7-H3, B7RP-2, 4lg-B7-H3) for WB, IHC-P, ELISA. Detects 57 kDa membrane protein involved in T-cell immune regulation.

[Learn More](#)

| Parameter        | Value                                                                                                                           |
|------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD276/B7-H3 Rabbit pAb                                                                                                          |
| Remarks/Alias    | B7H3; B7-H3; B7RP-2; 4lg-B7-H3; CD276/B7-H3                                                                                     |
| Species          | Human                                                                                                                           |
| GeneID (Human)   | 80381                                                                                                                           |
| GeneID           | 80381                                                                                                                           |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 29-245 of human CD276/B7-H3 (NP_079516.1). |
| Source           | Rabbit                                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                                           |
| Application      | WB, IHC-P, ELISA                                                                                                                |
| Cross Reactivity | Human                                                                                                                           |
| SWISS            | Q5ZPR3                                                                                                                          |
| Protein Weight   | 57kDa                                                                                                                           |
| Shipping         | Ice bag                                                                                                                         |