


Anti-Cd28 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P10747

Item Number: CM0002743



Introduction

CD28 Rabbit mAb validated for Western blot, flow cytometry, and ELISA. Recognizes human, mouse, and rat CD28; ~25 kDa. Essential tool for T-cell costimulation and immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	CD28 Rabbit mAb
Remarks/Alias	Tp44; CD28
Species	Human
GeneID (Human)	940
GeneID	940
Immunogen	Recombinant protein of human CD28 Full length P10747 .
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10747
Protein Weight	25kDa
Shipping	Ice bag