

Anti-Cd298/Atp1B3 Rabbit Polymab® For Wb, If/Icc, Ihc-P, Elisa - P54709

Item Number: CM0002410

Introduction

Rabbit monoclonal antibody against CD298/ATP1B3, a sodium/potassium-transporting ATPase subunit beta-3. Validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, and rat samples. Ideal for membrane transport research.

[Learn More](#)

Parameter	Value
Product Name	CD298/ATP1B3 Rabbit PolymAb®
Remarks/Alias	CD298; ATPB-3; ATP1B3
Species	Human
GeneID (Human)	483
GeneID	483
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 180-279 of human ATP1B3 (P54709)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54709
Protein Weight	22kDa/32kDa
Shipping	Ice bag