

Anti-Cd2Ap Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y5K6

Item Number: CM0027818



Introduction

CD2AP Rabbit pAb is a polyclonal antibody raised against recombinant human CD2AP, suitable for WB and ELISA, and cross-reactive with mouse and rat. Designed for research on cytoskeleton and cell junctions.

[Learn More](#)

Parameter	Value
Product Name	CD2AP Rabbit pAb
Remarks/Alias	CMS; CD2AP
Species	Human
GeneID (Human)	23607
GeneID	23607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-639 of human CD2AP (NP_036252.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5K6
Protein Weight	71kDa
Shipping	Ice bag