

Anti-Cd2Bp2 Polyclonal Antibody For Wb And Elisa - 095400

Item Number: CM0017508



Introduction

Rabbit polyclonal antibody targeting human CD2BP2 (U5-52K), validated for Western blot and ELISA. Involved in pre-mRNA splicing, suitable for studies of spliceosome assembly and U5 snRNP function.

[Learn More](#)

Parameter	Value
Product Name	CD2BP2 Rabbit pAb
Remarks/Alias	LIN1; Snu40; FWP010; U5-52K; PPP1R59; CD2BP2
Species	Human
GeneID (Human)	10421
GeneID	10421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-341 of human CD2BP2 (NP_006101.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95400
Protein Weight	38kDa
Shipping	Ice bag