

Anti-Cd3 Epsilon Monoclonal Antibody For Wb, If-P, Ihc-P, Fc, Elisa - P07766

Item Number: CM0006537



Introduction

Rabbit monoclonal antibody targeting human CD3 epsilon chain (T3E/TCRE), validated for WB, IF-P, IHC-P, FC, and ELISA. Ideal for T-cell activation and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	CD3 Rabbit mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag