

Anti-Cd3 Rabbit Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0009746



Introduction

-conjugated rabbit monoclonal antibody targeting the CD3 epsilon chain (23 kDa) of human, cynomolgus, and rhesus T cells. Designed for flow cytometry, it enables immunophenotyping and T-cell activation studies in basic research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus, Rhesus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag