


Anti-Cd30 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P28908

Item Number: CM0002424



Introduction

High-quality rabbit monoclonal antibody against human CD30 (TNFRSF8). Validated for WB, IHC-P, IF/ICC, FC, and ELISA. Ideal for lymphoma research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD30 Rabbit mAb
Remarks/Alias	CD30; Ki-1; D1S166E
Species	Human
GeneID (Human)	943
GeneID	943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-379 of human CD30 (NP_001234.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P28908
Protein Weight	14 kDa/50 kDa/63 kDa
Shipping	Ice bag