


Anti-Cd30 Rabbit Polyclonal Antibody For Wb, Elisa - P28908

Item Number: CM0026903



Introduction

CD30 Rabbit Polyclonal Antibody (P28908) validated for WB and ELISA in human. Synthetic peptide immunogen from amino acids 500-595. Ideal for lymphoblast activation and growth regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CD30 Rabbit pAb
Remarks/Alias	CD30; Ki-1; D1S166E
Species	Human
GeneID (Human)	943
GeneID	943
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-595 of human CD30 (NP_001234.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P28908
Protein Weight	64kDa
Shipping	Ice bag