

Anti-Cd300A Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q9Ugn4

Item Number: CM0008587



Introduction

Rabbit monoclonal antibody against human CD300A (IRp60). Suitable for Western blot, flow cytometry, and ELISA. Recognizes an inhibitory receptor on NK cells, mast cells, and other immune cells, with multiple isoforms detected.

[Learn More](#)

Parameter	Value
Product Name	CD300a Rabbit mAb
Remarks/Alias	IRC1, IRC2, CLM-8, IRp60, IGSF12, CMRF35H, CMRF-35H, CMRF35-H, CMRF35H9, CMRF35-H9, IRC1/IRC2, CMRF-35-H9
Species	Human
GeneID	11314
Immunogen	Recombinant protein corresponding to amino acids 18-178 of human CD300A (NP_009192.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	Q9UGN4
Protein Weight	33kDa/20kDa/12kDa/16.5kDa
Shipping	Ice bag