


Anti-Cd300A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ugn4

Item Number: CM0017941



Introduction

Rabbit polyclonal antibody against human CD300A, validated for WB and ELISA. Targets a 33 kDa inhibitory receptor on NK cells, T cells, and mast cells. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD300A Rabbit pAb
Remarks/Alias	IRC1; IRC2; CLM-8; IRp60; IGSF12; CMRF35H; CMRF-35H; CMRF35-H; CMRF35H9; CMRF35-H9; IRC1/IRC2; CMRF-35-H9; CD300A
Species	Human
Gene ID (Human)	11314
Gene ID	11314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-180 of human CD300A (NP_009192.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UGN4
Protein Weight	33kDa
Shipping	Ice bag