


Anti-Cd300A/C Monoclonal Antibody () For Flow Cytometry - Q9Ugn4/Q08708

Item Number: CM0008811



Introduction

-conjugated rabbit monoclonal antibody targeting human CD300a/c for flow cytometry. Suitable for immune cell profiling and functional studies of this inhibitory receptor family. Validated in FC application.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD300 a/c mAb
Remarks/Alias	IRC1; IRC2; CLM-8; IRp60; IGSF12; CMRF35H; CMRF-35H; CMRF35-H; CMRF35H9; CMRF35-H9; IRC1/IRC2; CMRF-35-H9
Species	Human
Gene ID (Human)	11314/10871
Gene ID	11314/10871
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-178 of human CD300 a (NP_009192.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UGN4/Q08708
Protein Weight	33kDa
Shipping	Ice bag