


Anti-Cd300C Polyclonal Antibody For Wb, Elisa - Q08708

Item Number: CM0012957



Introduction

Rabbit polyclonal antibody targeting human CD300C (CMRF35-like molecule 6). Validated for WB and ELISA applications, with cross-reactivity to mouse and rat. Ideal for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD300C Rabbit pAb
Remarks/Alias	LIR; CLM-6; CMRF35; IGSF16; CMRF-35; CMRF35A; CMRF-35A; CMRF35A1; CMRF35-A1; CD300C
Species	Human
Gene ID (Human)	10871
Gene ID	10871
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08708
Protein Weight	25kDa
Shipping	Ice bag