


Anti-Cd300C Rabbit Monoclonal Antibody For If/Icc, Fc, Elisa - Q08708

Item Number: CM0004126



Introduction

Human CD300C rabbit monoclonal antibody validated for IF/ICC, flow cytometry, and ELISA. Recombinant immunogen spanning the extracellular domain (aa 21-183). Detects CD300C, a leukocyte immune receptor expressed on monocytes, neutrophils, T and B cells. Ideal for immunology research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	CD300C Rabbit mAb
Remarks/Alias	LIR; CLM-6; CMRF35; IGSF16; CMRF-35; CMRF35A; CMRF-35A; CMRF35A1; CMRF35-A1; CD300C
Species	Human
GeneID (Human)	10871
GeneID	10871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-183 of human CD300C (NP_006669.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	Q08708
Protein Weight	25kDa
Shipping	Ice bag