

Anti-Cd300Lb Polyclonal Antibody For Wb, Elisa - A8K4G0

Item Number: CM0020451

Introduction

Recombinant immunogen-based CD300LB rabbit polyclonal antibody for Western blot and ELISA applications. Recognizes human, mouse, and rat CD300LB, an activating immune receptor on myeloid cells. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD300LB Rabbit pAb
Remarks/Alias	CLM7; CLM-7; IREM3; TREM5; CD300b; IREM-3; TREM-5; CMRF35-A2; CD300LB
Species	Human
GeneID (Human)	124599
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human CD300LB (NP_777552.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A8K4G0
Protein Weight	23kDa
Shipping	Ice bag