


Anti-Cd301/Clec10A Polyclonal Antibody For Wb, Elisa - Q8Iun9

Item Number: CM0016849



Introduction

Rabbit polyclonal antibody targeting CD301/CLEC10A, a C-type lectin receptor that recognizes GalNAc-terminated glycans. Validated for WB and ELISA, cross-reacts with Human and Mouse, supporting immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CD301/CLEC10A Rabbit pAb
Remarks/Alias	HML; MGL; HML2; CD301; CLECSF13; CLECSF14; CD301/CLEC10A
Species	Human
GeneID (Human)	10462
GeneID	10462
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-180 of human CD301/CLEC10A (NP_878910.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IUN9
Protein Weight	35kDa
Shipping	Ice bag