


Anti-Cd304/Neuropilin-1 Rabbit Polymab® For Wb, If/Icc, Ihc-P, Elisa - P97333

Item Number: CM0012433



Introduction

Rabbit monoclonal Neuropilin-1 antibody (PolymAb®) validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat NRP1. Ideal for studying angiogenesis, neuronal development, and cardiovascular biology.

[Learn More](#)

Parameter	Value
Product Name	CD304/Neuropilin-1 Rabbit PolymAb®
Remarks/Alias	Nrp; NP-1; Npn1; NPN-1; C530029I03
Species	Mouse
GeneID (Human)	8829
GeneID	18186
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P97333
Protein Weight	103kDa
Shipping	Ice bag