

Anti-Cd305/Lair1 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q6Gtx8

Item Number: CM0004765



Introduction

Rabbit monoclonal antibody against human leukocyte-associated immunoglobulin-like receptor 1 (LAIR1/CD305). Validated for WB, FC, ELISA. Ideal for immune receptor signaling and NK/T/B cell function studies.

[Learn More](#)

Parameter	Value
Product Name	CD305/LAIR1 Rabbit mAb
Remarks/Alias	LAIR1; CD305; LAIR-1; leukocyte-associated immunoglobulin-like receptor 1; CD305/LAIR1
Species	Human
GeneID (human)	3903
GeneID	3903
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-163 of human CD305/LAIR1(NP_002278.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	Q6GTX8
Protein Weight	23kDa/29kDa/31kDa
Shipping	Ice bag