

Anti-Cd307B Rabbit Monoclonal Antibody For Wb, If/Icc, Fc, Elisa - Q96La5

Item Number: CM0011662

Introduction

High-quality rabbit monoclonal antibody against human CD307b (FCRL2), ideal for WB, IF/ICC, FC, and ELISA applications in B-cell development research.

[Learn More](#)

Parameter	Value
Product Name	CD307b Rabbit mAb
Remarks/Alias	FCRH2; IFGP4; IRTA4; SPAP1; CD307b; SPAP1A; SPAP1B; SPAP1C
Species	Human
GeneID (Human)	79368
GeneID	79368
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 20-395 of human CD307b (NP_110391.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	Q96LA5
Protein Weight	56kDa
Shipping	Ice bag