

Anti-Cd31 Monoclonal Antibody () For Flow Cytometry - Q3Swt0

Item Number: CM0004212

Introduction

Rabbit monoclonal anti-rat CD31 antibody, conjugate, validated for flow cytometry. Detects platelet endothelial cell adhesion molecule (PECAM1) for leukocyte migration and endothelial adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Rat CD31 mAb
Remarks/Alias	CD31; Pecam; CD31/PECAM1
Species	Rat
Gene ID (Human)	5175
Gene ID	29583
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 18-585 of rat CD31 (NP_113779.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	Q3SWT0
Protein Weight	76kDa
Shipping	Ice bag