


Anti-Cd317 Rabbit Polyclonal Antibody For Wb, Elisa - Q8R2Q8

Item Number: CM0021904



Introduction

Rabbit polyclonal antibody against mouse CD317/Bst2 (bone marrow stromal antigen 2) for WB and ELISA applications. Suitable for studying innate antiviral immunity and viral restriction mechanisms.

[Learn More](#)

Parameter	Value
Product Name	CD317 Rabbit pAb
Remarks/Alias	GREG; Bst-2; CD317; DAMP-1; 2310015I10Rik
Species	Mouse
GeneID	69550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-152 of mouse CD317 (NP_932763.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8R2Q8
Protein Weight	19kDa
Shipping	Ice bag