

Anti-Cd31/Pecam1 Monoclonal Antibody (Apc/Cyanine7) For Flow Cytometry - Q08481

Item Number: CM0006103

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse/rat CD31/PECAM1. Validated for flow cytometry. Ideal for studying endothelial cell adhesion, leukocyte transmigration, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse
Gene ID (Human)	5175
Gene ID	18613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	69kDa/77kDa/80kDa/81kDa
Shipping	Ice bag