

Anti-Cd31/Pecam1 Monoclonal Antibody (Biotin) For Fc - Q08481

Item Number: CM0008253

Introduction

Biotin-conjugated rabbit monoclonal antibody against mouse/rat CD31/PECAM1 for flow cytometry. Validated in mouse and rat. Useful for studying endothelial cell adhesion, leukocyte transmigration, and platelet biology.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse
GenelD (Human)	5175
GenelD	18613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	69kDa/77kDa/80kDa/81kDa
Shipping	Ice bag