

Anti-Cd31/Pecam1 Mouse Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P16284

Item Number: CM0004743



Introduction

Monoclonal antibody targeting human CD31/PECAM1, validated for WB, IF-P, IHC-P, and ELISA. Shows cross-reactivity with mouse. Suitable for endothelial cell adhesion, leukocyte transmigration, and vascular biology studies.

[Learn More](#)

Parameter	Value
Product Name	CD31/PECAM1 Mouse mAb
Remarks/Alias	CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM
Species	Human
Gene ID (Human)	5175
Gene ID	5175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 630-738 of CD31/PECAM1 (NP_000433.4).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16284
Protein Weight	83kDa
Shipping	Ice bag