

Anti-Cd31/Pecam1 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q08481

Item Number: CM0026968

Introduction

Rabbit polyclonal antibody targeting CD31/PECAM1, validated for WB, IF-P, IHC-P, and ELISA. Detects mouse, human, and rat PECAM1, a critical endothelial adhesion molecule for leukocyte migration and angiogenesis. Ideal for vascular inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	CD31/PECAM1 Rabbit pAb
Remarks/Alias	Cd31; Pecam; PECAM-1; CD31/PECAM1
Species	Mouse
Gene ID (Human)	5175
Gene ID (Mouse)	18613
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 630-738 of CD31/PECAM1 (NP_000433.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08481
Protein Weight	81kDa
Shipping	Ice bag