

Anti-Cd31/Pecam1 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P16284

Item Number: CM0002930



Introduction

Rabbit monoclonal antibody against human CD31/PECAM1, validated for WB, IF, IHC, and ELISA. Detects human, mouse, and rat CD31. Suitable for studies on endothelial cell adhesion, leukocyte transmigration, and phagocytosis.

[Learn More](#)

Parameter	Value
Product Name	CD31/PECAM1 Rabbit mAb
Remarks/Alias	CD31; PECA1; GPIIA ⁺ ; PECAM-1; endoCAM; CD31/EndoCAM; CD31/PECAM1
Species	Human
Gene ID (Human)	5175
Gene ID	5175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 630-738 of CD31/PECAM1 (NP_000433.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16284
Protein Weight	83kDa
Shipping	Ice bag