

Anti-Cd31/Pecam1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P - P16284

Item Number: CM0014836



Introduction

High-quality rabbit polyclonal antibody targeting human CD31/PECAM1, suitable for Western blot, paraffin immunofluorescence, and immunohistochemistry. Reacts with human and mouse, supporting studies in leukocyte migration, endothelial junctions, and inflammatory diseases.

[Learn More](#)

Parameter	Value
Product Name	CD31/PECAM1 Rabbit pAb
Remarks/Alias	CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM
Species	Human
Gene ID	5175
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P
Cross Reactivity	Human, Mouse
SWISS	P16284
Protein Weight	83 kDa
Shipping	Ice bag