


Anti-Cd320 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npf0

Item Number: CM0025332



Introduction

High-quality rabbit polyclonal antibody against CD320 (TCbIR/8D6 antigen), validated for WB and ELISA in human and mouse samples. Suitable for studying B cell differentiation, proliferation, and cobalamin uptake.

[Learn More](#)

Parameter	Specification
Product Name	CD320 Rabbit pAb
Remarks/Alias	8D6; 8D6A; TCBLR; TCN2R; TCII-R; sCD320; CD320
Species	Human
Gene ID	51293
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-229 of human CD320 (NP_057663.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NPF0
Protein Weight	29kDa
Shipping	Ice bag