


## Anti-Cd324/E-Cadherin Polyclonal Antibody For Wb, Elisa - P09803

Item Number: CM0023352



### Introduction

Rabbit polyclonal antibody targeting mouse CD324/E-Cadherin, validated for WB and ELISA. Cross-reacts with human. Suitable for studying cell adhesion, epithelial biology, and cancer. 98 kDa protein. - P09803

[Learn More](#)

| Parameter        | Value                                                                                                                          |
|------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD324/E-Cadherin Rabbit pAb                                                                                                    |
| Remarks/Alias    | CDH1; Arc-1; CD324; CDHE; ECAD; LCAM; UVO; AA960649; ARC-1; E-cad; Ecad; L-CAM; Um; cadherin-1; CD324/E-Cadherin               |
| Species          | Mouse                                                                                                                          |
| GeneID (Human)   | 999                                                                                                                            |
| GeneID           | 12550                                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 157-369 of mouse CD324/E-Cadherin (NP_033994.1). |
| Source           | Rabbit                                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                                          |
| Application      | WB, ELISA                                                                                                                      |
| Cross Reactivity | Human, Mouse                                                                                                                   |
| SWISS            | P09803                                                                                                                         |
| Protein Weight   | 98kDa                                                                                                                          |
| Shipping         | Ice bag                                                                                                                        |