


Anti-Cd32A + Cd32B + Cd32C Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P12318 / P31994 / P31995

Item Number: CM0003805



Introduction

High-quality rabbit monoclonal antibody against human CD32A, CD32B, and CD32C. Validated for WB, IHC-P, and ELISA applications. Suitable for research on Fcγ receptors and immune response.

[Learn More](#)

Parameter	Value
Product Name	CD32A + CD32B + CD32C Rabbit mAb
Remarks/Alias	CD32; FCG2; FcGR; CD32A; CDw32; FCGR2; IGFR2; FCGR2A1; FcγR1a; CD32A + CD32B + CD32C
Species	Human
Gene ID	2212/2213/9103
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P12318/P31994/P31995
Protein Weight	35kDa/34kDa
Shipping	Ice bag