

Anti-Cd33 Rabbit Monoclonal Antibody For Wb, Elisa - P20138

Item Number: CM0007357



Introduction

Rabbit monoclonal antibody targeting human CD33 (SIGLEC3), validated for Western blot and ELISA applications. Recombinant immunogen (aa 18-259). Ideal for myeloid cell and microglia research in immunology and neuroinflammation.

[Learn More](#)

Parameter	Value
Product Name	CD33 Rabbit mAb
Remarks/Alias	p67; SIGLEC3; SIGLEC-3; CD33
Species	Human
GeneID (Human)	945
GeneID	945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-259 of human CD33 (NP_001763.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P20138
Protein Weight	40kDa
Shipping	Ice bag