

Anti-Cd33 Rabbit Polyclonal Antibody For Wb, Elisa - P20138

Item Number: CM0031059



Introduction

CD33 Rabbit Polyclonal Antibody specifically detects human CD33 (Siglec-3), a myeloid cell surface antigen. Validated for Western blot and ELISA applications, it is suitable for studies in immunology and neurobiology. Swiss-Prot: P20138.

[Learn More](#)

Parameter	Value
Product Name	CD33 Rabbit pAb
Remarks/Alias	p67; SIGLEC3; SIGLEC-3
Species	Human
GeneID (human)	945
GeneID	945
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 265-364 of human CD33 (NP_001763.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P20138
Protein Weight	40kDa
Shipping	Ice bag