

Anti-Cd335/Nkp46 Monoclonal Antibody For If/Icc, Fc, Elisa - O76036

Item Number: CM0007480



Introduction

Human CD335/NKp46 rabbit monoclonal antibody validated for IF/ICC, FC, ELISA; targets natural cytotoxicity triggering receptor 1 on NK cells; ideal for immune-oncology research.

[Learn More](#)

Parameter	Value
Product Name	CD335/NKp46 Rabbit mAb
Remarks/Alias	NCR1; CD335; LY94; NK-p46; NKP46; natural cytotoxicity triggering receptor 1
Species	Human
GeneID (Human)	9437
GeneID	9437
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-254 of human CD335/NKp46 (NP_004820.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	O76036
Protein Weight	21kDa/22kDa/23kDa/32kDa/34kDa
Shipping	Ice bag