

Anti-Cd335/Nkp46 Rabbit Polyclonal Antibody For Wb And Elisa - Q8C567

Item Number: CM0019538

Introduction

High-quality rabbit polyclonal antibody against mouse CD335 (NKp46, Ly94) validated for Western blot and ELISA. Ideal for NK cell cytotoxicity research. Swiss-Prot: Q8C567.

[Learn More](#)

Parameter	Value
Product Name	CD335 Rabbit pAb
Remarks/Alias	Ly94, NKp46, CD335
Species	Mouse
GeneID	17086 (mouse)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-255 of mouse CD335 (NP_034876.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8C567
Protein Weight	37kDa
Shipping	Ice bag