

Anti-Cd34 Monoclonal Antibody For If/Icc, Fc, Elisa - Q64314

Item Number: CM0006528

Introduction

Rabbit monoclonal antibody against mouse CD34, ideal for IF/ICC, flow cytometry, and ELISA applications in hematopoietic stem cell research. This antibody assists in studying cell adhesion and early hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	CD34 Rabbit mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Mouse
GeneID (human)	947
GeneID	12490
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-287 of mouse CD34 (NP_598415.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC, ELISA
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	35kDa/41kDa
Shipping	Ice bag