

Anti-Cd34 Rabbit Monoclonal Antibody For Wb, Ihc-P, Mihc, Elisa - P28906

Item Number: CM0003540

Introduction

Recombinant rabbit monoclonal antibody against human CD34. Validated for WB, IHC-P, mIHC, ELISA; cross-reacts with human, mouse, rat. Suitable for hematopoietic stem and endothelial cell research.

[Learn More](#)

Parameter	Value
Product Name	CD34 Rabbit mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Human
GeneID (Human)	947
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 286-385 of human CD34 (P28906)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28906
Protein Weight	41 kDa
Shipping	Ice bag