

Anti-Cd34 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28906

Item Number: CM0017782

Introduction

High-quality CD34 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human CD34, a hematopoietic progenitor cell marker involved in cell adhesion. Ideal for stem cell and endothelial research.

[Learn More](#)

Parameter	Value
Product Name	CD34 Rabbit pAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Human
GeneID (Human)	947
GeneID	947
Immunogen	Synthetic peptide, a synthetic peptide corresponding to a sequence within amino acids 32-290 of human CD34 (NP_001020280.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P28906
Protein Weight	41 kDa
Shipping	Ice bag