

Anti-Cd34 Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P28906

Item Number: CM0010291

Introduction

Anti-CD34 Rabbit PolymAb® monoclonal antibody, validated for WB, IHC-P, IF/ICC, and ELISA in human samples. Suitable for hematopoietic stem cell and endothelial research, detecting 41/35 kDa isoforms.

[Learn More](#)

Parameter	Value
Product Name	CD34 Rabbit PolymAb®
Remarks/Alias	CD34
Species	Human
GeneID (Human)	947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P28906
Protein Weight	41kDa/35kDa
Shipping	Ice bag