

Anti-Cd35 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P17927

Item Number: CM0007907



Introduction

High-quality rabbit monoclonal antibody against human CD35/CR1 (P17927) for complement research. Validated for WB, IHC-P, IF/ICC, and ELISA; targets immune adherence receptor in innate immunity.

[Learn More](#)

Parameter	Value
Product Name	CD35 Rabbit mAb
Remarks/Alias	KN; C3BR; C4BR; CD35
Species	Human
Gene ID (Human)	1378
Gene ID	1378
Immunogen	Recombinant protein corresponding to amino acids 901-1095 of human CD35 (NP_000564.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P17927
Protein Weight	224kDa
Shipping	Ice bag