

Anti-Cd351/Fcamr Monoclonal Antibody For Elisa, If/Icc, Fc - Q8Wwv6

Item Number: CM0009493

Introduction

High affinity rabbit monoclonal antibody targeting human CD351/FCAMR (Fc alpha/mu receptor). Validated for ELISA, IF/ICC, and flow cytometry, suitable for immune response and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	CD351/FCAMR Rabbit mAb
Remarks/Alias	CD351; FCA/MR; FKSG87; Fcalpha/muR
Species	Human
GenelD (human)	83953
GenelD	83953
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 17-450 of human CD351/FCAMR (Q8WWV6-1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	ELISA, IF/ICC, FC
Cross Reactivity	Human
SWISS	Q8WWV6
Protein Weight	24kDa/46kDa/57kDa
Shipping	Ice bag