

Anti-Cd352/Slamf6 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q9Et39

Item Number: CM0009073

Introduction

Recombinant protein-generated rabbit monoclonal antibody targeting mouse CD352/SLAMF6 (Ly108). Designed for WB, IF-P, ELISA. Ideal for analyzing SLAM family signaling, lymphocyte activation, and immune regulation in mouse hematopoietic cells.

[Learn More](#)

Parameter	Value
Product Name	CD352/SLAMF6 Rabbit mAb
Remarks/Alias	KAL1; NTBA; KAL1b; Ly108; NTB-A; SF2000
Species	Mouse
GeneID (Human)	114836
GeneID	30925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-235 of mouse CD352/SLAMF6 (NP_001334115.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q9ET39
Protein Weight	39kDa
Shipping	Ice bag