

Anti-Cd352/Slamf6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Du3

Item Number: CM0007633



Introduction

Rabbit monoclonal antibody against human CD352/SLAMF6 (SLAM family member 6), validated for WB, IHC-P, and ELISA. Detects a 24/37 kDa band. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD352/SLAMF6 Rabbit mAb
Remarks/Alias	KALI; NTBA; CD352; KALib; Ly108; NTB-A; SF2000; CD352/SLAMF6
Species	Human
GeneID	114836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-215 of human CD352/SLAMF6 (NP_001171643.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q96DU3
Protein Weight	24 kDa/37 kDa
Shipping	Ice bag