

Anti-Cd353/Slamf8 Monoclonal Antibody For Flow Cytometry - Q9P0V8

Item Number: CM0010428



Introduction

-conjugated rabbit monoclonal antibody against human SLAMF8 (CD353, BLAME) for flow cytometry. Explore B-cell receptor signaling and lineage commitment with reliable IVD raw materials.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD353/SLAMF8 mAb
Remarks/Alias	BLAME; CD353; SBB142
Species	Human
Gene ID	56833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-123 of human CD353/SLAMF8 (NP_064510.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9P0V8
Protein Weight	32kDa
Shipping	Ice bag